

Mitochondria segmentation in EM data with Napari-Empanada

Marie Held, PhD
Image Analyst

Segmentation of Electron microscopy images

- Electron microscopy images stained with heavy metals
- Grey scale images containing dozens of organelles & subcellular structures

YES

Highly integrated
understanding re. biological
function

BUT

Major image analysis
challenge!

- Specific identification of different organelle types in electron microscopy is based on some combination of the electron density, contrast, shape, and relative location of each structure

Empanada

1. **Purpose:** provide essential nutrients, energy for bodily functions, and support for growth, repair, and overall health
2. **Key Feature:** type of baked or fried turnover consisting of pastry and filling, common in Spain, other Southern European countries, North African countries, Latin American countries, and the Philippines.
3. **Dataset:**



Napari-Empanada

- 1. Purpose:** Empanada-napari plugin democratizes deep learning image segmentation for electron microscopy (EM) researchers
- 2. Key Feature:** Comes pre-loaded with **MitoNet**, a pre-trained specialist machine learning model for mitochondria instance segmentation in EM images.
- 3. Dataset:**
 - Pre-trained on ~1.5 million unlabeled cellular EM images.
 - Fine-tuned on ~135,000 labeled mitochondrial instances.

Napari-Empanada Background

1. **Documentation:**

https://empanada.readthedocs.io/en/latest/getting_started/install.html

2. **MitoNet publication:**

<https://www.sciencedirect.com/science/article/pii/S240547122200494X>

Installation

1. Install Anaconda Navigator
2. Install Napari

```
conda update -n base conda  
conda create -n empanada-ws-env python=3.9  
conda activate empanada-ws-env  
conda install -c conda-forge napari pyqt
```

3. Install Empanada-napari plugin

Note: do NOT install plugin through Napari GUI

```
pip install empanada-napari
```

Usage – Case 1

1. Demo Dataset: <https://mitoem.grand-challenge.org/MitoEM/>
 - ~ 26 GB → tiny crop
2. 3D Inference:
 - Segmentation Confidence Thr: 0.50
 - Centre Confidence Thr: 0.10
 - Centres Min Distance: 3

Usage – Case 2

1. Demo Dataset 2: Cell Gatan 3View Serial Block Face Cell

2. Load data

3. Scale data as voxels are anisotropic:

```
viewer.layers[Image_name'].scale = [5, 1, 1]
```

4. 3D Inference:

- Segmentation Confidence Thr: 0.35
- Centre Confidence Thr: 0.75
- Centres Min Distance: 5

Usage – Case 2

- Segmentation is ok but room for improvement
→ Pick fine tuning patches
 - # Patches for annotation: 1
 - Patch size (px) = 224
 - Delete initial dataset and labels
 - Use 2D inference for a “head start”
 - Get initial segmentation for “clean up”
 - Tick [Output to layer]
 - On third slice of flip book
 - Run 2D inference → segmentations output to labels layer

Usage – Case 2


- Paint and erase labels on third slice only
 - Duplicate labels layer (if required)
 - Make edits on original labels layer
 - Select label using dropper [4]
 - Erase [1], draw [2], fill [3] as required
 - See [documentation](#) for merging, splitting of labels
- Save training patches and labels

Usage – Case 2

- Fine tune model – ideally on GPU leveraging hardware
 - Training directory: annotated images
 - Finetunable layers:
 - **None**: Select this option if the model did fairly well during the initial inference on your data.
 - **Stages 1 – 4**: Select between these options depending on how well the model did on the initial inference.
 - **All**: Select this option if the model did not perform well on your data. This option will take more time but could offer better results depending on your specific task.

Summary

- Plugin designed to simplify the application of deep learning models for electron microscopy (EM) image segmentation.
- suited for researchers aiming to perform instance segmentation of mitochondria, using pre-trained model MitoNet
- Ideal for large-scale vEM data sets
- Leverage pre-trained deep learning models without coding experience
- Fine tune pre-trained models to custom model



MHELDB@LIVERPOOL.AC.UK

[HTTP://CCI.LIV.AC.UK](http://CCI.LIV.AC.UK)

